

Product Datasheet

Tenalisib R Enantiomer (orb1180502)

Description	Tenalisib R Enantiomer
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{23}H_{18}FN_5O_2$
Purity	> 98%
MW	415.42
CAS Number	[1639417-54-1]
SMILES	<chem>O=C1C(C2=CC=CC(F)=C2)=C([C@H](NC3=C4N=CNC4=NC=N3)CC)OC5=CC=CC=C15</chem>
Expiration Date	12 months from date of receipt.