

Product Datasheet

3,4-Benzocoumarin (orb1180516)

Description	3,4-Benzocoumarin
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	C ₁₃ H ₈ O ₂
Purity	> 98%
MW	196.2
CAS Number	[2005-10-9]
SMILES	<chem>O=C1C2=CC=CC=C2C3=CC=CC=C3O1</chem>
Expiration Date	12 months from date of receipt.