

Product Datasheet

DPC-423 (orb1182670)

Description	DPC-423
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{25}H_{21}ClF_4N_4O_3S$
Purity	> 98%
MW	568.97
CAS Number	[292135-59-2]
SMILES	<chem>CS(=O)(=O)C1=CC=CC=C1C2=CC(=C(C=C2)NC(=O)C3=CC(=NN3C4=CC=CC(=C4)CN)C(F)(F)F)</chem>
Expiration Date	12 months from date of receipt.