

Product Datasheet

A-841720 (orb1182733)

Description	A-841720
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{17}H_{21}N_5OS$
Purity	> 98%
MW	343.45
CAS Number	[869802-58-4]
SMILES	<chem>CN(C)C1=C2C3=C(C(=O)N(C=N3)N4CCCCC4)SC2=NC=C1</chem>
Expiration Date	12 months from date of receipt.