

Product Datasheet

trans-Zeatin (orb1183368)

Description	trans-Zeatin
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{10}H_{13}N_5O$
Purity	> 98%
MW	219.25
CAS Number	[1637-39-4]
SMILES	<chem>OC/C(C)=C/CNC1=NC=NC2=C1N=CN2</chem>
Expiration Date	12 months from date of receipt.