

Product Datasheet

Strontium ranelate (Protelos) (orb1183452)

Description	Strontium ranelate (Protelos)
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	C ₁₂ H ₆ N ₂ O ₈ S.2Sr
Purity	> 98%
MW	513.49
CAS Number	[135459-87-9]
SMILES	C(C1=C(SC(=C1C#N)N(CC(=O)[O-])CC(=O)[O-])C(=O)[O-])C(=O)[O-].[Sr+2].[Sr+2]
Expiration Date	12 months from date of receipt.