

Product Datasheet

gamma-secretase modulator 2 (orb1183507)

Description	gamma-secretase modulator 2
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{25}H_{22}F_4N_6O_2$
Purity	> 98%
MW	514.5
CAS Number	[1093978-89-2]
SMILES	<chem>CC1=CN(C=N1)C2=C(C=C(C=C2)C3=CN(N=N3)C4CCC5=C(C=CC=C5F)N(C4=O)CC(F)(F)F)OC</chem>
Expiration Date	12 months from date of receipt.