

Product Datasheet

MK-8245 (orb1183514)

Description	MK-8245
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{17}H_{16}BrFN_6O_4$
Purity	> 98%
MW	467.3
CAS Number	[1030612-90-8, 1030612-87-3]
SMILES	<chem>C1CN(CCC1OC2=C(C=CC(=C2)F)Br)C3=NOC(=C3)C4=NN(N=N4)CC(=O)O</chem>
Expiration Date	12 months from date of receipt.