

Product Datasheet

EMD638683 S-Form (orb1183515)

Description	EMD638683 S-Form
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{18}H_{18}F_2N_2O_4$
Purity	> 98%
MW	364.34
CAS Number	[1184940-46-2]
SMILES	<chem>CC1=C(CC)C(C(NNC([C@@H](O)C2=CC(F)=CC(F)=C2)=O)=O)=CC=C1O</chem>
Expiration Date	12 months from date of receipt.