

Product Datasheet

YM 750 (orb1183517)

Description	YM 750
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{31}H_{36}N_2O$
Purity	> 98%
MW	452.63
CAS Number	[138046-43-2]
SMILES	<chem>O=C(NC1=C(C)C=C(C)C=C1C)N(CC2=CC=C3C(CC4=CC=CC=C34)=C2)C5CCCCCCC5</chem>
Expiration Date	12 months from date of receipt.