

Product Datasheet

Dipraglurant (orb1183519)

Description	Dipraglurant
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	C ₁₆ H ₁₂ FN ₃
Purity	> 98%
MW	265.29
CAS Number	[872363-17-2]
SMILES	<chem>FC1=CN2C(C=C1)=NC(CCC#CC3=NC=CC=C3)=C2</chem>
Expiration Date	12 months from date of receipt.