

Product Datasheet

Calcitetrolo (orb1183523)

Description	Calcitetrolo
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	C ₂₇ H ₄₄ O ₄
Purity	> 98%
MW	432.64
CAS Number	[72203-93-1]
SMILES	<chem>C[C@H](CCC(C(C)(C)O)O)[C@H]1CC[C@@H]2[C@@]1(CCC/C2=CC=C/C[C@H](C[C@@H](C3=C)O)O)C</chem>
Expiration Date	12 months from date of receipt.