

Product Datasheet

CB 300919 (orb1183525)

Description	CB 300919
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{32}H_{34}ClN_7O_2$
Purity	> 98%
MW	584.11
CAS Number	[289715-28-2]
SMILES	<chem>CN1CCN(CC1)CC2=NC3=CC(=C(C=C3C(=O)N2C)CN(CC#C)C4=CC=C(C=C4)C(=O)NCC5=CN=CC=C5)Cl</chem>
Expiration Date	12 months from date of receipt.