

Product Datasheet

Ciluprevir (BILN 2061) (orb1183530)

Description	Ciluprevir (BILN 2061)
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	C ₄₀ H ₅₀ N ₆ O ₈ S
Purity	> 98%
MW	774.9
CAS Number	[300832-84-2]
SMILES	<chem>CC(C)NC1=NC(=CS1)C2=NC3=C(C=CC(=C3)OC)C(=C2)O[C@@H]4C[C@@H]5C(=O)N[C@@]6(C[C@H]6/C=CCCCC[C@@H](C(=O)N5C4)NC(=O)OC7CCCC7)C(=O)O</chem>
Expiration Date	12 months from date of receipt.