

Product Datasheet

BMS 433796 (orb1183532)

Description	BMS 433796
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{19}H_{16}F_2N_4O_4$
Purity	> 98%
MW	402.4
CAS Number	[935525-13-6]
SMILES	<chem>CN1C(=O)[C@H](C2=CC=CC=C2C=N1)NC(=O)NC(=O)[C@H](C3=CC(=CC(=C3)F)F)O</chem>
Expiration Date	12 months from date of receipt.