

Product Datasheet

P005672 HCl (Sarecycline HCl) (orb1183542)

Description	P005672 HCl (Sarecycline HCl)
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{24}H_{30}ClN_3O_8$
Purity	> 98%
MW	524
CAS Number	[1035979-44-2]
SMILES	<chem>CN(C)[C@H]1[C@@H]2C[C@@H]3CC4=C(C=CC(=C4C(=C3C(=O)[C@@]2(C(=C(C1=O)C(=O)N)O)O)O)O)CN(C)OC.Cl</chem>
Expiration Date	12 months from date of receipt.