

Product Datasheet

AEG 3482 (orb1183543)

Description	AEG 3482
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{10}H_8N_4O_2S_2$
Purity	> 98%
MW	280.33
CAS Number	[63735-71-7]
SMILES	<chem>NS(C1=NN(C(S1)=N2)C=C2C3=CC=CC=C3)(=O)=O</chem>
Expiration Date	12 months from date of receipt.