

Product Datasheet

Cenisertib (orb1183545)

Description	Cenisertib
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{24}H_{30}FN_7O$
Purity	> 98%
MW	451.54
CAS Number	[871357-89-0]
SMILES	<chem>[H][C@@]1([C@H]([C@H]2C(N)=O)NC3=NC(NC4=CC=C(C(C)=C4)N5CCN(CC5)C)=NC=C3F)C[C@]2(C=C1)[H]</chem>
Expiration Date	12 months from date of receipt.