

Product Datasheet

PF-03084014 (orb1183546)

Description	PF-03084014
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{27}H_{41}F_2N_5O$
Purity	> 98%
MW	489.64
CAS Number	[865773-15-5, 1290543-63-3]
SMILES	<chem>CCC[C@@H](C(=O)NC1=CN(C=N1)C(C)(C)CNCC(C)(C)C)N[C@H]2CCCC3=CC(=CC(=C3C2)F)F</chem>
Expiration Date	12 months from date of receipt.