

Product Datasheet

XL647 (Tesevatinib) (orb1183583)

Description	XL647 (Tesevatinib)
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{24}H_{25}Cl_2FN_4O_2$
Purity	> 98%
MW	491.4
CAS Number	[874286-84-7]
SMILES	<chem>CN1CCCC(CCC1)COC2=C(C=C3C(=C2)N=CN=C3NC4=CC(=C(C(=C4)Cl)Cl)F)OC</chem>
Expiration Date	12 months from date of receipt.