

Product Datasheet

MK-4827 (Niraparib) (orb1183589)

Description	MK-4827 (Niraparib)
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{19}H_{20}N_4O$
Purity	> 98%
MW	320.4
CAS Number	[1038915-60-4]
SMILES	<chem>C1C[C@H](CNC1)C2=CC=C(C=C2)N3C=C4C=CC=C(C4=N3)C(=O)N</chem>
Expiration Date	12 months from date of receipt.