

Product Datasheet

TAK-875 (Fasiglifam) (orb1183590)

Description	TAK-875 (Fasiglifam)
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	C ₂₉ H ₃₂ O ₇ S
Purity	> 98%
MW	524.6
CAS Number	[1000413-72-8]
SMILES	<chem>CC1=CC(=CC(=C1C2=CC(=CC=C2)COC3=CC4=C(C=C3)[C@@H](CO4)CC(=O)O)C)OCCCS(=O)(=O)C</chem>
Expiration Date	12 months from date of receipt.