

## Product Datasheet

### MAP3K12 antibody (orb1183774)

|                            |                                       |
|----------------------------|---------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to MAP3K12 |
| <b>Species/Host</b>        | Rabbit                                |
| <b>Reactivity</b>          | Human, Mouse, Rat                     |
| <b>Conjugation</b>         | Unconjugated                          |
| <b>Tested Applications</b> | ELISA, ICC, WB                        |
| <b>Target</b>              | MAP3K12                               |
| <b>Preservatives</b>       | PBS                                   |
| <b>Form/Appearance</b>     | Liquid                                |
| <b>Concentration</b>       | 1mg/ml                                |
| <b>Note</b>                | For research use only                 |
| <b>Clonality</b>           | Polyclonal                            |
| <b>Expiration Date</b>     | 12 months from date of receipt.       |