

## Product Datasheet

### GCLM antibody (orb1183823)

|                            |                                    |
|----------------------------|------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to GCLM |
| <b>Species/Host</b>        | Rabbit                             |
| <b>Reactivity</b>          | Human, Mouse, Rat                  |
| <b>Conjugation</b>         | Unconjugated                       |
| <b>Tested Applications</b> | ELISA, ICC, IHC                    |
| <b>Target</b>              | GCLM                               |
| <b>Preservatives</b>       | PBS                                |
| <b>Form/Appearance</b>     | Liquid                             |
| <b>Concentration</b>       | 1mg/ml                             |
| <b>Note</b>                | For research use only              |
| <b>Clonality</b>           | Polyclonal                         |
| <b>Expiration Date</b>     | 12 months from date of receipt.    |