


Product Datasheet

CD68 antibody (orb1183943)

Description	CD68 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human CD68
Target	CD68
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	37 kDa
Uniprot ID	P34810
Dilution Range	IHC-P=1:20-1000, Flow-Cyt=1ug/Test, ICC=1:50, IF=1:50
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
human tonsil
using CD6...