


Product Datasheet

CD74 antibody (orb1183946)

Description	CD74 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant mouse CD74 protein (56-279/279 aa)
Target	CD74
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	34 kDa
Uniprot ID	P04233
Dilution Range	IHC-P=1:400-800
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
mouse
thymus using
CD7...