

Product Datasheet

CD178 antibody (FITC) (orb1184086)

Description	Mouse monoclonal antibody CD178 conjugated to FITC
Species/Host	Mouse
Conjugation	FITC
Tested Applications	FC
Target	FASLG
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Form/Appearance	Liquid
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Isotype	Mouse IgG1, k
Clonality	Monoclonal
MW	31 kDa
Uniprot ID	AAH17502
Dilution Range	Flow-Cyt=
Expiration Date	12 months from date of receipt.