

Product Datasheet

CD42b antibody (PE/Cy5) (orb1184138)

Description	Mouse monoclonal antibody CD42b conjugated to PE/Cy5
Species/Host	Mouse
Conjugation	PE/Cy5
Tested Applications	FC
Target	GP1BA
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Form/Appearance	Liquid
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Isotype	Mouse IgG1, k
Clonality	Monoclonal
MW	67 kDa
Uniprot ID	P07359
Dilution Range	Flow-Cyt=
Expiration Date	12 months from date of receipt.