

Product Datasheet

CD123 antibody (PE) (orb1184171)

Description	Mouse monoclonal antibody CD123 conjugated to PE
Species/Host	Mouse
Conjugation	PE
Tested Applications	FC
Target	IL3RA
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Form/Appearance	Liquid
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Isotype	Mouse IgG1, k
Clonality	Monoclonal
MW	40 kDa
Uniprot ID	P26951
Dilution Range	Flow-Cyt=
Expiration Date	12 months from date of receipt.