

Product Datasheet

CD45.2 antibody (PE) (orb1184177)

Description	Mouse monoclonal antibody CD45.2 conjugated to PE
Species/Host	Mouse
Conjugation	PE
Tested Applications	FC
Target	PTPRC
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Form/Appearance	Liquid
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Isotype	Mouse IgG2a, k
Clonality	Monoclonal
MW	180-240 kDa
Uniprot ID	P08575
Dilution Range	Flow-Cyt=
Expiration Date	12 months from date of receipt.