

Product Datasheet

CD20 antibody (PE) (orb1184203)

Description	Rat monoclonal antibody CD20 conjugated to PE
Species/Host	Rat
Conjugation	PE
Tested Applications	FC
Target	MS4A1
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Form/Appearance	Liquid
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Isotype	Rat IgG2a, k
Clonality	Monoclonal
MW	33 kDa
Uniprot ID	P11836
Dilution Range	Flow-Cyt=
Expiration Date	12 months from date of receipt.