

Product Datasheet

CD45.2/Pacific Blue antibody (orb1184225)

Description	Mouse monoclonal antibody to CD45.2/Pacific Blue
Species/Host	Mouse
Conjugation	Pacific Blue
Tested Applications	FC
Target	PTPRC
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Form/Appearance	Liquid
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Isotype	Mouse IgG2a, k
Clonality	Monoclonal
MW	180-240 kDa
Uniprot ID	P08575
Expiration Date	12 months from date of receipt.