

Product Datasheet

CD49e antibody (PE) (orb1184233)

Description	Mouse monoclonal antibody CD49e conjugated to PE
Species/Host	Mouse
Conjugation	PE
Tested Applications	FC
Target	ITGA5
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Form/Appearance	Liquid
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Isotype	Rat IgG2a, k
Clonality	Monoclonal
MW	110 kDa
Uniprot ID	P08648
Dilution Range	Flow-Cyt=
Expiration Date	12 months from date of receipt.