

Product Datasheet

CD80 antibody (FITC) (orb1184241)

Description	Rat monoclonal antibody CD80 conjugated to FITC
Species/Host	Rat
Conjugation	FITC
Tested Applications	FC
Target	CD80
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Form/Appearance	Liquid
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Isotype	Rat IgG2a, k
Clonality	Monoclonal
MW	32 kDa
Uniprot ID	P33681
Expiration Date	12 months from date of receipt.