

Product Datasheet

Ki-67 antibody (FITC) (orb1184302)

Description	Mouse monoclonal antibody Ki-67 conjugated to FITC
Species/Host	Mouse
Conjugation	FITC
Tested Applications	FC
Target	MKI67
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Form/Appearance	Liquid
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	358 kDa
Uniprot ID	P46013
Dilution Range	Flow-Cyt=
Expiration Date	12 months from date of receipt.