

Product Datasheet

CD11b/Pacific Blue antibody (orb1184317)

Description	Rat monoclonal antibody to CD11b/Pacific Blue
Species/Host	Rat
Conjugation	Pacific Blue
Tested Applications	FC
Target	ITGAM
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Form/Appearance	Liquid
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Isotype	Rat IgG2b, κ
Clonality	Monoclonal
Uniprot ID	P11215
Expiration Date	12 months from date of receipt.