

Product Datasheet

CD62L (L-Selectin) antibody (FITC) (orb1184359)

Description	Rat monoclonal antibody CD62L (L-Selectin) conjugated to FITC
Species/Host	Rat
Conjugation	FITC
Tested Applications	FC
Target	SELL
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Form/Appearance	Liquid
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Isotype	Rat IgG2a, k
Clonality	Monoclonal
MW	88 kDa
Uniprot ID	P14151
Dilution Range	Flow-Cyt=
Expiration Date	12 months from date of receipt.