

Product Datasheet

CD107a antibody (PE/Cy5) (orb1184439)

Description	Mouse monoclonal antibody CD107a conjugated to PE/Cy5
Species/Host	Mouse
Conjugation	PE/Cy5
Tested Applications	FC
Target	LAMP1
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Form/Appearance	Liquid
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Isotype	Mouse IgG1, k
Clonality	Monoclonal
MW	42 kDa
Uniprot ID	P11279
Expiration Date	12 months from date of receipt.