

Product Datasheet

CD107a/Pacific Blue antibody (orb1184445)

Description	Mouse monoclonal antibody to CD107a/Pacific Blue
Species/Host	Mouse
Conjugation	Pacific Blue
Tested Applications	FC
Target	LAMP1
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Form/Appearance	Liquid
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Isotype	Mouse IgG1, k
Clonality	Monoclonal
MW	42 kDa
Uniprot ID	P11279
Expiration Date	12 months from date of receipt.