

Product Datasheet

CD86 antibody (FITC) (orb1184485)

Description	Rat monoclonal antibody CD86 conjugated to FITC
Species/Host	Rat
Reactivity	Human
Conjugation	FITC
Tested Applications	FC
Target	CD86
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Form/Appearance	Liquid
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Isotype	Rat IgG2a, κ
Clonality	Monoclonal
MW	31 kDa
Uniprot ID	O35531
Dilution Range	Flow-Cyt=10ul/Test
Expiration Date	12 months from date of receipt.