

Product Datasheet

CD29 antibody (FITC) (orb1184629)

Description	Mouse monoclonal antibody CD29 conjugated to FITC
Species/Host	Mouse
Conjugation	FITC
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from human CD29
Target	ITGB1
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Isotype	Mouse IgG1, k
Clonality	Monoclonal
MW	88 kDa
Uniprot ID	P05556
Dilution Range	Flow-Cyt=
Expiration Date	12 months from date of receipt.