



## Product Datasheet

### NDUFV2 Antibody (orb1184710)


<b>Description</b>	NDUFV2 Antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, FC, IF, IHC, WB
<b>Immunogen</b>	E.coli-derived human NDUFV2 recombinant protein (Position: M1-L249).
<b>Form/Appearance</b>	Lyophilized
<b>Concentration</b>	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
<b>Storage</b>	At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.
<b>Note</b>	For research use only
<b>Application notes</b>	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
<b>Isotype</b>	Rabbit IgG
<b>Clonality</b>	Polyclonal
<b>MW</b>	27 kDa
<b>Uniprot ID</b>	<b>P19404</b>
<b>Dilution Range</b>	Western blot, 0.1-0.25 µg/ml/ml, Human, Mouse, Rat Immunohistochemistry(Paraffin-embedded Section), 2-5 µg/ml/ml, Human Flow Cytometry, 1-3 µg/ml/1x10 <sup>6</sup> cells, Human Direct ELISA, 0.1-0.5 µg/ml/ml, Human
<b>Expiration Date</b>	12 months from date of receipt.



IHC analysis of NDUFV2 using anti-NDUFV2...



IHC analysis of NDUFV2 using anti-NDUFV2...



IHC analysis of NDUFV2 using anti-NDUFV2...