





Product Datasheet

IMPDH2 Antibody (orb1184715)

Description	IMPDH2 Antibody	
Species/Host	Rabbit	 IHC analysis of IMPDH2 using anti-IMPDH2...
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	FC, ICC, IF, IHC, WB	
Immunogen	A synthetic peptide corresponding to a sequence of human IMPDH2 (VVLSPKDRVRDVFEAKARH).	 IHC analysis of IMPDH2 using anti-IMPDH2...
Form/Appearance	Lyophilized	
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.	
Storage	At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.	
Note	For research use only	 IHC analysis of IMPDH2 using anti-IMPDH2...
Application notes	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.	
Isotype	Rabbit IgG	
Clonality	Polyclonal	
MW	56 kDa	 IHC analysis of IMPDH2 using anti-IMPDH2...
Uniprot ID	P12268	
Dilution Range	Western blot, 0.25-0.5 µg/ml/ml, Human Immunohistochemistry(Paraffin-embedded Section), 2-5 µg/ml/ml, Human Immunocytochemistry/Immunofluorescence, 5 µg/ml, Human Flow Cytometry, 1-3 µg/ml/1x10 ⁶ cells, Human	
Expiration Date	12 months from date of receipt.	