

Product Datasheet

CPSF6 Antibody (monoclonal, 3F11E1) (orb1184752)

Description

CPSF6 Antibody (monoclonal, 3F11E1)

Species/Host

Mouse

Reactivity

Human, Monkey, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

FC, ICC, IF, IHC, WB

Immunogen

E.coli-derived human CPSF6 recombinant protein (Position: R50-Q176).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Isotype

Mouse IgG1

Clonality

Monoclonal

MW

68 kDa

Uniprot ID


[Q16630](#)

Dilution Range


Western blot, 0.25-0.5 µg/ml, Human, Mouse, Rat, Monkey
Immunohistochemistry(Paraffin-embedded Section), 2-5 µg/ml, Human, Mouse, Rat
Immunocytochemistry/Immunofluorescence, 5 µg/ml, Human Flow Cytometry, 1-3 µg/1x10⁶ cells, Human

Expiration Date


12 months from date of receipt.



IHC analysis of CPSF6 using anti-CPSF6 a...



IHC analysis of CPSF6 using anti-CPSF6 a...



IHC analysis of CPSF6 using anti-CPSF6 a...