

Product Datasheet

NKCC1/SLC12A2 Antibody (monoclonal, 6G7D2) (orb1184753)

Description

NKCC1/SLC12A2 Antibody (monoclonal, 6G7D2)

Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IF, IHC, WB

Applications
Immunogen

E.coli-derived human NKCC1/SLC12A2 recombinant protein (Position: K889-K943).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Isotype

Mouse IgG2a

Clonality

Monoclonal

MW


200 kDa

Uniprot ID
P55011
Dilution Range


Western blot, 0.25-0.5 µg/ml, Human Immunohistochemistry(Paraffin-embedded Section), 2-5 µg/ml, Human, Mouse, Rat Immunofluorescence, 5 µg/ml, Human, Mouse, Rat

Expiration Date


12 months from date of receipt.



IHC analysis of NKCC1/SLC12A2 using anti...



IHC analysis of NKCC1/SLC12A2 using anti...



IHC analysis of NKCC1/SLC12A2 using anti...