

Product Datasheet

Human KLRG1 Protein (orb1184764)

Description	Human KLRG1, His Tag is expressed from human 293 cells (HEK293). It contains AA Leu 60 - Phe 195 (Ac...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	KLRG1
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	N-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the N-terminus. The protein has a calculated MW of 17.4 kDa. The protein migrates as 30-42 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_001316028.1
Purity	90%
Source	Human KLRG1, His Tag (orb1184764) is expressed from human 293 cells (HEK293). It contains AA Leu 60 - Phe 195 (Accession # Q96E93-1).
MW	17.4 kDa
Uniprot ID	Q96E93
NCBI	NP_001316028.1
Expiration Date	6 months from date of receipt.