

## Product Datasheet

### Human IL-28 R alpha / IFN-lambda R1 Protein (orb1184778)

**Description**

Human IL-28 R alpha, His Tag is

(HEK293). It contains AA Arg 21 - Ala...

**Reactivity**

Human

**Endotoxins**

1.0 EU per µg

**Conjugation**

Unconjugated

**Target**

IL-28 R alpha

**Preservatives**

PBS, pH7.4

**Form/Appearance**

Powder

**Storage**

-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

**Tag**

C-10×His

**Note**

For research use only

**Application notes**

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 25.5 kDa. The protein migrates as 40-45 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

**Purity**

95%

**Source**

Human IL-28 R alpha, His Tag (orb1184778) is expressed from human 293 cells (HEK293). It contains AA Arg 21 - Ala 228 (Accession # Q8IU57-1).

**MW**

25.5 kDa

**Uniprot ID**

**Q8IU57**

**NCBI**

**NP\_734464.1**

**Expiration Date**

6 months from date of receipt.