

Product Datasheet

**Biotinylated Human CCL1 / I-309 / C-C motif
chemokine 1 Protein (orb1184779)**

Description

Biotinylated Human CCL1, His, is

(HEK293). It contains AA Lys 24 - Ly...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Biotin

Target

CCL1

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Tag

C-10×His & C-Avi

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus, followed by an Avi tag (). The protein has a calculated MW of 12.1 kDa. The protein migrates as 15 kDa and 17-19 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

90%

Source

Biotinylated Human CCL1, His, (orb1184779) is expressed from human 293 cells (HEK293). It contains AA Lys 24 - Lys 96 (Accession # P22362-1).

MW

12.1 kDa

Uniprot ID

P22362

NCBI

NP_002972.1

Expiration Date

6 months from date of receipt.