

Product Datasheet

Human NKp46 / NCR1 / CD335 Protein (orb1184795)

Description	Natural cytotoxicity triggering receptor 1 (NCR1) is also known as Natural killer cell p46-related p...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	Nkp46
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 28.9 kDa. The protein migrates as 35-45 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_004820.2(K82Q)
Purity	90%
Source	Human NKp46, His Tag (orb1184795) is expressed from human 293 cells (HEK293). It contains AA Gln 22 - Arg 258 (Accession # O76036-1).
MW	28.9 kDa
NCBI	K82, NP_004820.2
Expiration Date	6 months from date of receipt.