

Product Datasheet

CD303 antibody (APC) (orb1185304)

Description	Mouse monoclonal antibody against CD303 conjugated to APC
Species/Host	Mouse
Reactivity	Human
Conjugation	APC
Tested Applications	FC
Preservatives	PBS with 0.05% Proclin300, 1% rAlbumin
Storage	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Mouse IgG2a, κ
Uniprot ID	Q8WTT0
Entrez	170482
Expiration Date	12 months from date of receipt.